




BOD and CBOD

PolySeed® and PolySeed NX®

- *Easy-to-use capsules*
- *Convenient and cost-effective BOD analysis with consistent results*
- *Free of nitrifying microorganisms*
- *EPA recommended*

PolySeed® The “Original” BOD Seed

Have fresh seed on demand every day without handling wastewater. Each capsule of PolySeed contains 100 mg of specialized microbial cultures ideal for use in a broad range of areas involving degradation of industrial and municipal waste. Get convenient, fast, and economical BOD analysis with consistent results using InterLab’s PolySeed.



PolySeed NX® The “Original” CBOD Seed

Obtain accurate CBOD readings even in the presence of a healthy nitrifying population. A blend of specialized microbial cultures containing a chemical additive for the inhibition of nitrifying bacteria in wastewater samples. Designed to provide a standard seed source for the degradation of both industrial and municipal waste in CBOD testing.


Biomass Toxicity

PolyTox®

- *Determine toxicity to wastewater biomass*
- *Reliable results in 30 minutes*
- *Non-pathogenic*
- *Simple and economical*
- *EPA recommended*

PolyTox® Rapid Biological Toxicity Screening Test

Rapid test for measuring toxicity of pH, temperature, specific chemicals, individual influent streams or wastewater to biological systems. A blend of specialized microbial cultures in an easy-to-use kit provides a simple and economical test.



The Market Leader in Biological Wastewater Testing Products

For more information on BOD and toxicity testing, please visit www.polyseed.com. As the market leader in biological wastewater testing products, InterLab Supply’s consumer promise is to provide the highest quality product, exceptional technical support without exception and prompt, accurate delivery.

Product Characteristics:

PolySeed: Granular powder form
Buff to brown color
Musty odor

PolySeed NX: Granular powder form
Buff to brown color
Musty odor

PolyTox: Granular powder form
Buff to brown color
Musty odor

Application Instructions:

Please contact us directly or visit our website at www.polyseed.com for complete product application information and ordering options.

Packaging:

Polyseed: 50 capsule bottles
Polytox: Package of 20 vials

Polyseed NX: 50 capsule bottles

Safety:

InterLab products contain naturally occurring, nongenetically engineered bacteria which have been subjected to stringent antibiotic screening. Simple hygiene practices are recommended with use. Do not ingest this product. Direct contact with the skin should be avoided if irritation is noted. In the event of direct skin or eye contact, flush the affected area with water. If irritation persists, contact your physician.

Storage:

Polyseed, Polyseed NX, Polytox: Store sealed in dry area at room temperature between 65°F and 80°F (18°-20°C).

Contact:

InterLab Supply®
Technical Service Department
4200 Research Forest Drive, Suite 150
The Woodlands, TX 77381
Phone: 888.MyPolySeed
email: polyseed@polyseed.com
www.polyseed.com