

PolySeed[®]

Convenient, fast and cost-effective BOD₅ analysis with consistent results. Each capsule contains a minimum of 100 mg of specialized microbial cultures ideal for use in a broad range of areas involving degradation of industrial and municipal waste. Formulated according to *The Standard Methods for Examination of Water and Wastewater*.

- EPA accepted
- Non-pathogenic
- Free of nitrifying microorganisms
- To use, mix contents of one (1) capsule with 0.5L of APHA standard dilution water at 20°C (68°F), stir for at least 60 minutes
- Mixture yields seed inoculum for up to 200 BOD₅ tests
- Shelf life: 2 years
- Also available as **PolySeed NX**, which contains an approved chemical additive for the inhibition of ammonia nitrogen (NH₃-N) in wastewater samples for CBOD₅ analysis